



Repair Services - EDL Chiropractic Products



As the manufacturer of Nervo-Scope® and Analagraph® Chiropractic instruments - we pride ourselves in our capability to repair our products efficiently and economically.

If the need for repair arises, you can forward your instrument to us for a free evaluation and repair quotation. Upon your approval, your instrument is repaired, refurbished and tested to insure performance of optimum functionality.

We continually strive for the fastest turnaround to limit our customer's inconvenience - as we have done for over 50 years! The staff at Electronic Development Labs, Inc. is happy to address your specific questions and concerns.

Chiropractic Products

Nervo-Scope®
Thermal Detection
Instruments

Analagraph®
Graph Recording
Instrument



Electronic Development Labs, Inc

244 Oakland Drive
Danville, VA 24540, USA
1-800-DIAL-EDL
Ph: 434-799-0807
Fax: 434-799-0847
website: <http://www.nervo-scope.com>



Electronic Development Labs, Inc

244 Oakland Drive
Danville, VA 24540, USA
1-800-DIAL-EDL
Ph: 434-799-0807
Fax: 434-799-0847
website: www.edl-inc.com
email: sales@edl-inc.com



Nervo-Scope® - Thermal Detection Instruments

Nervo-Scope® reliably detects and measures minute heat differences in the tissues along the spine. The instrument can be used comfortably for both pre and post adjustment assessment.

Nervo-Scope® is built to withstand the rigors of continual cleaning and disinfecting. A large meter scale with a fluorescent pointer makes even the smallest temperature change easy to read. The sensing detectors can be adjusted from 0.6" wide to 2.1" wide to accommodate infants and adults.

Each Nervo-Scope® includes a foam-padded carrying case and cleaning brush. All of the ETS series Nervo-Scopes are protected by a full two-year warranty against workmanship and defective components.



ETS-6

Nervo-Scope® ETS-6 is a stand alone version. The ETS-6 is for use in applications where chart recording is not necessary. The ETS-6 is not compatible with the Analagraph® Recorder.



ETS-7

Nervo-Scope® ETS-7 is for use as a stand alone unit or in applications where chart recording is necessary. The ETS-7 is compatible with the Analagraph® Recorder.



ETS-8

The exceptional stability of the Nervo-Scope® line has been improved by the introduction of the ETS-8. Instrument self heating, caused by the chiropractors hand or repeated use, will not dampen the response of the ETS-8.

Features

- **Meter** - Our taut-band meter construction assures continual durability and reliability.
- **Large Scale** - For fast readings with large pointer for easy tracking.
- **Comfort Glide Caps** - Self lubricated for smooth, easy glides.
- **Highly Sensitive Detectors** - Capable of picking up temperature changes as little as 0.02°F.
- **Adjustable Detector Width** - Capable of adjusting from 0.6" wide to 2.1" wide.
- **Function Reliability** - All components used are selected for their long-term reliable performance and stability.
- **Auto-Power On/Off** - One easy-touch button powers the instrument. After three minutes of non-use, the unit will automatically power off.
- **Multiple Sensitivity Settings** - The ETS-6 and ETS-7 features six sensitivity settings. The ETS-8 offers two sensitivity settings.
- **Durable Construction** - Easy-use design molded from mineral-filled nylon for durability, chemical resistance and non-slip surface.

About the ETS-8

Because of the highly developed circuitry used in the ETS-8, only two sensitivities are required. The ETS-8 is highly stable and virtually unaffected by the heat generated during use. Signal amplitude is fixed to each sensitivity and is not dampened by the heat generated by the chiropractor's hand.

The ETS-8 is compatible with the Analagraph recording unit and is constructed of the highly reliable components used in the original ETS-6 and ETS-7 Nervo-Scopes. The unit features an all black case to distinguish it from the other models.

A great deal of engineering effort was put into the ETS-8 to achieve this improvement. It is an improvement that will benefit both the chiropractor and their patient.

Analagraph® - Graph Recorder

The Analagraph® Recorder, when used in conjunction with the Nervo-Scope® ETS-7 or ETS-8, provides a permanent record to log and chart the amount of existing heat differences indicated from your Nervo-Scope® readings.


The Analagraph® Recorder is a valuable and efficient tool in documenting patient assessments and providing a pattern reference for pre and post check of patients. This information indicates any change of reading, thus establishing a pattern which may assist in determining the frequency of adjustments.

Easy to use and operate, the Analagraph® Recorder has been trusted for over 50 years in the chiropractic community.

CHR1

(Shown with ETS-7 sold separately)

Power Input:	120Vac @ 60Hz, 1A
Overload Protection:	3 AG, 250V slo-blow glass body fuse (rear panel)
Ambient Operating Range:	32° to 140°F
Relative Humidity (max.):	10% to 90% non-condensing
Recorder Resolution:	0.02°F
Power Consumption:	Recorder - 5 watts @ 0.3 amps on recording mode. 1/2 watt on standby
Printing Method:	Direct thermal
Print Resolution:	200 x 800 DPI
Graph Speed Rate:	Low: 6 inches per minute (2.5 mm sec.) High: 9 inches per minute (3.75 mm sec.)
Graph Paper Type:	Thermal activation
Graph Paper Size:	2 inch wide roll x 100 feet long
Graph Paper Thickness:	0.0025"
Graph Grid Size:	0.195" squares x 8 horizontally x continuous vertically
Switching:Power:	On/Off rocker switch on rear panel
Graph On/Off:	Starts/stops thermal graphing
Time Marker:	Automatic horizontal dashed lines every 3 vertical grid squares
Manual Marker:	Allows manual placement of horizontal dashed lines at any vertical interval
High/Low Speed:	Sets graph paper speed rate
Audible Sound:	Beep occurs for "graph on", "time marker" & "manual marker" intervals
Green Light indicator:	For "graph on" (active)
Red Light Indicator:	For "graph off" (standby)
ETS-7/ETS-8 Input Socket:	Accepts coil cord, linking Nervo-Scope to Analagraph
Coil Cord:	6 ft. (expandable to 12 ft.) 6 wire pickup
Case:	High impact black ABS with handle - powder coated front/rear panels
Size:	8.5" width x 6.4" height x 9.2" deep - not including handle
Weight:	5 lbs.



**Replacement Recording Paper Rolls
#CHR-AGP-CASE**
Case / 100 Rolls
Single Rolls
(10 Roll Minimum)

Sporicidin® Disinfectant Products are ideal for industrial & pharmaceutical clean rooms, hospitals, medical offices, dental offices, and laboratories.



Sporicidin® Disinfectant Solutions

SSS-1284C	Sterilizing Solution 1 gal. size
RE-1284C	Disinfectant Solution and Presoak - 1 gal. size
PS-2212C	Disinfectant Pump Spray - 22 oz. size
CFS-0481	Disinfectant Pump Spray - 1 oz. (Qty. 36/box)

Sporicidin® Disinfectant Towelettes

ST50-600	Disinfectant Individually Wrapped Towelettes - 5" x 8", Qty. 50/box, Citrus Scent
CAN-18012C	Disinfectant Towelettes - 5" x 6-3/4", Qty. 160/canister, Citrus Scent
JUM-8506C	Jumbo size extra strength - 9.5" x 12", 85 towelettes/canister

Sporicidin® Antimicrobial Lotion Soap

CALS-2416	Antimicrobial Lotion Soap - 16 oz. pump
ALS-1282	Antimicrobial Lotion Soap - 1 gal. size, pump refill
WIPE-15144	Antimicrobial Handwipes - Qty. 15 Towelettes/packet (resealable), 12 Packets/box, Citrus Scent